



SHRADDHA
ANALYTICAL SERVICES

Serving The Industry Since Last 22 Years

Thermogravimetry Analysis (TGA)

Sample Requirement :-

1. Solid Sample - 1 gm .
2. Liquid Sample - 1 ml.



Applications :-

1. For measurement of Weight loss with respect to temp.
2. For measurement of Volatile Matter [Water, Solvent etc.]
3. Filler Content of Polymer Samples.

contactus - DR. KETAN MERCHANT [M.Sc., Ph.D - ICT]

Send your Enquiries



saslab21@gmail.com
info@saslabindia.com



www.saslabindia.com



+91 8850413449
+91 022 2500 0186
+91 022 2500 0386
+91 022 4215 0502


Sample: 2202821-1-CHO 098-033
Size: 17.7530 mg
Method: Ramp
Comment: AIR

DSC-TGA

File: C:\...\2202821-1-CHO 098-033.001

Run Date: 19-May-2022 09:24

Instrument: SDT Q600 V20.9 Build 20



Sample: 2202821-2-CHO 098-037
Size: 5.7580 mg
Method: Ramp
Comment: AIR

DSC-TGA

File: C:\...2202821-2-CHO 098-037.001

Run Date: 19-May-2022 11:36

Instrument: SDT Q600 V20.9 Build 20

